

Dip & Read Tests

For Science

Rapid Tests for Science

- Quick & easy
- High informative value
- Accurate results



MACHERY-NAGEL

www.mn-net.com

MN







Since 1911

MN Dip & Read Tests for Science

Dip & Read tests allow the easy detection of important parameters directly at the point of interest. They provide reliable results within a few seconds and allow a quick decision on further activities. MN test strips are the ideal products for everyday lab use and are successfully used in science labs around the globe.



Product information

Test	Product	Application	Gradation	REF
	pH-Fix 0-14	<ul style="list-style-type: none"> Always ready for use Patented technology prevents bleeding of the indicator Brilliant colour charts for exact results 	pH 0 · 1 · 2 · 3 · 4 · 5 · 6 · 7 · 8 · 9 · 10 · 11 · 12 · 13 · 14	921 10
	pH-Fix 0-14 PT	<ul style="list-style-type: none"> pH-Fix tests also available in PlopTop tubes Virtually unbreakable pack Easy to open and close 	pH 0 · 1 · 2 · 3 · 4 · 5 · 6 · 7 · 8 · 9 · 10 · 11 · 12 · 13 · 14	921 11
	QUANTOFIX® Peroxide 100	<ul style="list-style-type: none"> Peroxide containing solvents may explode when heated Test allows to check for peroxide in a few seconds Easy Dip & Read procedure for best usability 	0 · 1 · 3 · 10 · 30 · 100 mg/l H ₂ O ₂	913 12
	Potassium iodide starch paper	<ul style="list-style-type: none"> Detects strong oxidisers (i.e. free chlorine) Easy control of diazotization reactions (excess of nitrite) 	Detection limit 1 mg/l NO ₂ ⁻ or free Cl ₂	907 54
	Fluoride Test Paper	<ul style="list-style-type: none"> Detects dangerous hydrofluoric acid without additional reagents Easy control of puddles with unknown liquids Provides immediate results 	Detection limit 20 mg/l F ⁻	907 50
	Watesmo	<ul style="list-style-type: none"> Test paper for detection of water Reaction colour irreversible Easy control whether solvents are anhydrous 	Detection of water droplets	906 09

Contact

Your local distributor



Image Credits © Leah-Anne Thompson; Tom; NIDerLander - Fotolia.com

MACHERY-NAGEL

MACHERY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · D-52355 Düren · Germany
 Germany and international: Tel.: +49 (0) 24 21 96 90 Fax: +49 (0) 24 21 96 91 99 e-mail: sales-de@mn-net.com
 Switzerland: Tel.: +41 (0) 62 388 55 00 Fax: +41 (0) 62 388 55 05 e-mail: sales-ch@mn-net.com
 France: Tel.: +33 (0) 3 88 68 22 68 Fax: +33 (0) 3 88 51 76 88 e-mail: sales-fr@mn-net.com
 USA: Tel.: +1 484 821 0984 Fax: +1 484 821 1272 e-mail: sales-us@mn-net.com

